

TUMOUR MARKER CONTROL - LEVEL 2 (TMR CONTROL 2)

CAT. NO. TU5002

LOT NO. 310TU

SIZE: 3 x 2 ml

EXPIRY: 2022-06-28

GTIN: 05055273207828

INTENDED USE

This product is intended for *in vitro* diagnostic use, in the quality control of diagnostic assays on clinical chemistry and immunoassay systems. The Tumour Marker Controls are for the control of accuracy and reproducibility.

DEVICE DESCRIPTION

The Tumour Marker Controls are supplied at 2 levels, level 2 and 3. Target values and ranges are supplied for tumour markers, as listed in the value tables for both levels.

SAFETY PRECAUTIONS AND WARNINGS

For *in vitro* diagnostic use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.

Human source material, from which this product has been derived, has been tested at donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B Surface Antigen (HbsAg), and Hepatitis C Virus (HCV) antibody and found to be NON-REACTIVE. FDA approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, this material and all patient samples should be handled as though capable of transmitting infectious diseases and disposed of accordingly.

Health and Safety Data Sheets are available on request.

STORAGE AND STABILITY

OPENED: Store refrigerated (+2°C to +8°C). Once reconstituted, Tumour Marker Controls are stable for 14 days when stored tightly capped at +2°C to +8°C in the absence of contamination, with the following exceptions: Thyroglobulin and Calcitonin should be assayed immediately following reconstitution. No claim is made for the stability of CA 72-4, Calcitonin, Cyfra 21 and NSE. Only the required amount of product should be removed. After use, any residual product should NOT BE RETURNED to the original vial.

UNOPENED: Store refrigerated (+2°C to +8°C). Stable to expiration date printed on individual vials.

PREPARATION FOR USE

Open the vial carefully, avoiding any loss of the material and reconstitute with 2 ml of distilled water. Replace the rubber stopper, close the vial and leave to stand for 30 minutes before use. Ensure that all traces of dry material are dissolved by swirling gently.

MATERIALS PROVIDED

Tumour Marker Control - Level 2 3 x 2 ml

ASSIGNED VALUES

Each batch of Tumour Marker Control is submitted to a number of external laboratories and values are assigned from a consensus of results obtained by these laboratories. With each batch, a control range is provided for individual parameters and each parameter method.

EC REP

Randox Teoranta, Meenmore,
Dungloe, Donegal,
F94 TV06, Ireland

13 Nov '20 ne

TUMOUR MARKER CONTROL - LEVEL 2 (TMR CONTROL 2)

Cat. No. TU5002		Lot No. 310TU		Size: 3 x 2 ml		Expiry: 2022-06-28	
Analyte	unit	Target	Range		methods		
			low	high			
Alpha-fetoprotein	KIU/l = IU/ml	23.3	18.6	28.0	Roche Cobas Systems		
	ng/ml	28.2	22.6	33.8			
CA 15-3	U/ml	46.1	36.9	55.3	Roche Cobas Systems		
CA 19-9	U/ml	14.2	11.4	17.0	Roche Cobas Systems		
CA 72-4	U/ml	4.57	3.43	5.71	Roche Cobas Systems		
CA125	U/ml	67.3	53.8	80.8	Roche Cobas Systems		
Carcinoembryonic Antigen (CEA)	ng/ml = µg/l	6.14	4.91	7.37	Roche Cobas Systems		
Cyfra 21-1	ng/ml	5.00	3.75	6.25	Roche Cobas Systems		
Neuron Specific Enolase (NSE)	ng/ml	1.84	1.38	2.30	Roche Cobas Systems		
Thyroglobulin	ng/ml	11.6	8.70	14.5	Roche Cobas Systems		
Total beta hCG	mIU/ml = IU/l	6.04	4.83	7.25	Roche Cobas Systems		
	IU/ml	0.006	0.005	0.007			