

# TOTAL ANTIOXIDANT STATUS (TAS CONTROL) CONTROL SERUM

**CAT. NO.** NX 2331

**LOT NO.** 316NX

**SIZE:** 5 ml

**EXPIRY:** 2018-07

This Control Serum is human in origin.

## PREPARATION

Open the vial very carefully, avoiding any loss of material and reconstitute in exactly 5 ml of freshly distilled water at +20°C to +25°C. Replace the rubber stopper, close vial and leave to stand for 30 minutes before use.

## STABILITY

Stable for 2 days at +2°C to +8°C or 12 hours at +15°C to +25°C.

Target value (mmol/l)	Range (mmol/l)
1.72	1.38 - 2.06

## SAFETY

Human source material, from which this product has been derived, has been tested at donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B Surface Antigen (HbsAg), and Hepatitis C Virus (HCV) antibody and found to be NON-REACTIVE. FDA approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, this material and all patient samples should be handled as though capable of transmitting infectious diseases and disposed of accordingly.

07 Jan '15 ne

**THIS PAGE IS INTENTIONALLY BLANK**